实验室多功能提取浓缩回收机组 LAB MULTI-FUNCTION EXTRACTING TANK

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Pilotech experimental multi-function extraction thickening recycling set is a new developed product based on reference to domestic and foreign similar products, suitable for extraction of Chinese herbal medicine, natural product, health care product, biological pharmacy, cosmetics and food, together with supersonic extraction, ultrasonic extraction and dynamic extraction under normal pressure and temperature, extract vacuum relief concentration (water and organic solvent) as well as volatile aromatic oil recovery, with a concentration proportion of 1.1-1.38. This equipment can achieve dynamic and continuous ultrasonic extraction and vacuum relief concentration at ultra-low temperature. The equipment comprises of extraction tank, reflux condenser, reflux cooler, aromatic oil collector, extraction tank filter, tubular heater, evaporator, recovery condenser, recovery cooler, gathering tank, stainless steel frame, connection pipe valve and instrument, perfectly suitable for varied and small batch production and application in research and development department of universities, research institutes, enterprises and pharmaceutical factories.

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Equipment function. Static water extraction or alcohol extraction, hot reflux water extraction or alcohol extraction under normal pressure, external reflux concentration under normal pressure, concentration under normal pressure, vacuum relief concentration, organic solvent recovery (under normal pressure or reduced pressure) (online ultrasonic extraction function available)

设备特点 Features

- 1) 占地面积小, 仅需1~3M², 接上水和电就可使用, 最小投料量0.5KG, 最大投料量5KG。
Small area coverage of 2-3m². Works directly with water and power with a minimum inventory of 0.5KG and a maximum inventory of 5KG.
- 2) 提取系统配有回流冷凝器、回流冷却器、油水分离器, 可有效的回收物料中挥发性芳香油。
Extraction system is equipped with a reflux condenser, reflux cooler, oil-water separator that can recover aromatic oil in materials in an effective manner.
- 3) 夹层加热(夹层无压力), 加热均匀, 应用于有机溶剂提取浓缩更安全。
Circular heating of conduction oil (no pressure between lamination) and uniform heating. The heater is far away from the extraction tank which is more secure for organic solvent extraction and concentration.
- 4) 机组自带蒸发浓缩器, 浓缩效率更高, 浓缩比重可达1.1~1.38。
Units are equipped with outer circular evaporating concentrator with a more concentrating efficiency, concentration ratio reaching 1.1~1.38.
- 5) 可适应水提、醇提或溶剂提取, 是天然植物药提取方面的中试设备, 适用于科研机构、大专院校、医院、药厂、食品行业对中草药液、茶汁和果汁等物质的提取、浓缩以及有机溶剂的回收。
Can be applied to water, alcohol and solvent extraction and serve as pilot plant testing equipment for natural herb extraction. Suitable for extraction and concentration of herbal liquid and tea and fruit juice and recovery of organic solvents in research institutions, universities, hospitals, pharmaceutical factory and food industry.
- 6) 它集多功能提取、真空浓缩与有机溶剂回收为一体, 可进行热回流循环提取, 浓缩, 常规提取, 浓缩等操作, 简化了操作程序。
It combines multi-functional extraction, vacuum concentration and organic solvent recovery as a whole and simplify operations with functions of heat reflux circular extraction, concentration, conventional extraction, concentration, water sedimentation (alcohol sedimentation).
- 7) 机组与物料接触部分全部采用不锈钢制造, 能在负压, 常压状态下工作。Unit parts contacting materials are made of stainless steel that works in negative pressure, ordinary pressure and positive pressure conditions.
- 8) 低温浓缩, 保证产品质量, 动态循环蒸发, 提高了换热效率, 有效解决换热效率随浓度增高而降低问题, 节约能耗, 提高工效。
Low temperature concentration ensures product quality. Dynamic circular evaporation increases heat transfer efficiency, provides an effective solution to decreased heat transfer efficiency with increased concentration, saves energy and improves efficiency.
- 9) 药渣通过浸提篮提升出料, 极大的改善工作环境和提升工作效率。
Dregs leach through the digestion basket to greatly improve the working environment and efficiency.
- 10) 选配在线超声波发生器多功能提取罐特点
Features of online ultrasonic generator multi-function extraction tank

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High efficiency. On the basis of the original dynamic hot reflux extraction set, the newly developed superminiature ultrasonic dynamic extraction and condensing set optimizes its structure by combining energy-gathering ultrasonic generator with miniature vacuum relief condenser, which can complete ultrasonic dynamic extraction, abstraction, filter, relief concentration and essential oil condensation in one step, greatly saving raw materials and work time, with an improved work efficiency of 60% - 80% than general multi-function extraction machine.

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High raw material conversion rate. Unique physical function and cavitation of ultrasonic wave can promote plant cell tissue wall breaking or deformation and strengthen solute particle vibration, acceleration shock, sound pressure shear effect so that materials form local extreme high temperature and high pressure to facilitate traditional Chinese medicine effective component precipitation, with significantly improved extraction efficiency 45% higher than traditional processes.

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Ultrasonic extraction carries out under normal temperature, which effectively reduces heat-sensitive ingredient loss. In addition, ultrasonic extraction and condensing set adopts titanium alloy energy-gathering generator where big-size Chinese herbal medicine can fully contact with the ultrasonic probe to accelerate homogeneous precipitation of effective ingredients.

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The equipment has a compact structure and gives a full play to the ultrasonic energy-gathering generator. Ultrasonic waves directly act on materials, strengthens herbal extraction mechanical action, cavitation function, local high frequency vibration, high impact and high sound pressure shear effect to shorten the extraction time by 4/5, together with big volume of raw material processing.

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Ultrasonic extraction equipment and vacuum relief condenser can work independently and jointly.

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High raw material conversion rate. During concentration, evaporation chamber works at a lower temperature to retain heat-sensitive effective ingredients to the greatest degree. In this way, medicinal effect after concentration can be better retained.

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Wide adaptability. Chinese herbal medicine extraction is not restricted by ingredient polarity or molecular weight size, making it suitable for extraction of ingredients of most Chinese herbal medicine. This equipment is fully enclosed inner loop structure, equipped with oil-water separator, condenser and cooler, for repeated use on valuable solvents.

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